



Content

1. General information	3
2. Introduction	4
3. Display & operating buttons	5
4. Preparing before operating	6
5. Turn on & measure	7
6. Calibration	8
7. Changing scale & temperature unit	10
8. Turning off	11
9. Cleaning & maintenance	11
10. Disposal	11
11. Technical data	12
12. Error codes	12
13. Models and scales	13

⚠ Carefully read through the operating manual even if you have prior experience with KERN refractometers.

1. General information

1.1 Intended use

The refractometer is a measuring instrument for determining the refractive index of transparent substances in liquid or in some cases also in the solid state. It is used to observe the behaviour of light as it passes from a prism with known properties to the substance being tested. Use of the refractometer for other purposes is contrary to its intended use and may be hazardous. The manufacturer shall not be liable for any damages caused by improper use.

1.2 Warranty

The warranty shall be void in the event of:
 • Failure to observe the instructions in the operating manual
 • Use for purposes other than those described
 • Modifications or opening the device housing
 • Mechanical damage and/or damage resulting from media, liquids, natural wear and tear



This digital refractometer cannot measure any liquid that is highly corrosive to metal or glass. When measuring liquids that are corrosive to plastics or react chemically with plastics, be careful not to drop the measured liquid onto the shell. Otherwise it will corrode the shell.

2. Introduction


2.1 Description



2.2 Scope of delivery
 1x Storage box | 1x Digital refractometer | 1x Operating manual | 1x AAA Battery 1.5 V | 1x Pipette | 1x Screwdriver

3. Display & operating buttons



3.1 Description display & operating buttons



⚠ Note: Please replace the battery when the is displayed.

4. Preparing before operating


4.1 Install the battery



Put 1 piece of 1.5V battery into the cabin in the right way and recover the cabin again.
 Only for KERN service staff


5. Turn on & measure

5.1 Turn on




The multi-function display shows the current scale number. Press "READ" to turn on the refractometer. If there is no sample filled in, the display shows "- - -".

Note :
 1. When used outdoors, please avoid strong light so as not to affect the measurement accuracy.
 2. Please keep the instrument in a stable and still statement and position.



5.2 Measure

5.3 Average value measurement



Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.




Remaining measurements

6. Calibration


The refractometer can only be calibrated with distilled water. To do this, fill the sample tank with distilled water up to the mark and close the cover.

7. Changing scale & temperature unit


7.1 Changing scale



Press "CAL" for 2 seconds to enter calibration mode. Then press "CAL" again for 2-3 seconds until "CAL" flashes in the display.




While "CAL" is flashing in the display, press "CAL" again to start the calibration. When the calibration is finished, the display shows "End". After approx. 10 seconds, the device automatically returns to normal mode.




If the calibration was not completed successfully, an error code appears in the display. Here, for example, A01.


7.2 Changing temperature unit



Press "SCALE" to change into another scales and show the converted value.




To change the temperature unit, press "SCALE" for 2 seconds.



Further error codes can be found in the appendix.

8. Turning off

8.1 Turning off



If without any operations for 1 minute, the instrument would be automatically turned off.

9. Cleaning & maintenance

- To avoid damages to the prism and the sample tank, clean them with distilled water after each use.
- Dry it with a soft cloth afterwards.
- Do not use hard or abrasive objects for cleaning.
- Do not leave any residue in the sample tank.
- If the refractometer is not going to be used for a longer time, remove the battery and store it at a cool and dry place.

10. Disposal

The packaging consists of environmentally friendly materials which can be disposed of via local recycling facilities. The device and storage box should be disposed

11. Technical data

Scale + accuracy + resolution	Depents to the model
Temperature	0,0 – 40,0 °C / 32,0 – 104,0 °F
Automatic Temperature Compensation	Yes
Minimum sample volume	0.2 - 0.3 ml (Marking ring)
AUTO-OFF	60 seconds
Averaging measurement	15 measurements
Battery	1 x AAA 1.5 V
Lifetime of the battery	Approx. 10.000 measurements
Overall dimensions LxWxH	125x65x30 mm
Net weight	140 g (without battery)

12. Error codes

code	Instructions
A01	Beyond the scope of calibration temperature. (0.0°C-40.0°C)
A02	During calibration, no solution or solution wrong.
A03	This instrument has a hardware failure.

13. Models and scales

Model	Scale	No.	Range	Unit	Resolution	Accuracy
ORM 100M	Refractive Index	502	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Brix	501	0.0-99.9	%	0.1%	±0.2%
	Refractive Index	502	1.3330-1.5177	nD	0.0001nD	±0.0003nD
ORM 1RS	Refractive Index	502	0.0-99.9	%	0.1%	±0.2%
	Brix	501	0.0-99.9	%	0.1%	±0.2%
	Refractive Index	502	1.3330-1.5177	nD	0.0001nD	±0.0003nD
ORM 1SU	Fructose	501	0.0-88.9	%	0.1%	±0.2%
	Glucose	502	0.0-59.9	%	0.1%	±0.2%
	Brix	503	0.0-99.9	%	0.1%	±0.2%
ORM 2SU	Refractive Index	504	1.3330-1.5177	nD	0.0001nD	±0.0003nD
	Brix	501	0.0-18.5	%	0.1%	±0.2%
	Maltool	502	0.0-15.6	%	0.1%	±0.2%
ORM 1HD	Dehydran	503	0.0-19.6	%	0.1%	±0.2%
	Brix	504	0.0-50.0	%	0.1%	±0.2%
	Honey Water	501	0.0-38.0	%	0.1%	±0.2%
ORM 1HA	Honey Bouma	502	0.0-28.0	%	0.1%	±0.2%
	Brix	503	0.0-99.0	%	0.1%	±0.2%
	Refractive Index	504	1.3330-1.5177	nD	0.0001nD	±0.0003nD
ORM 1NA	Salinity NaCl (%)	501	0.0-28.0	%	0.1%	±0.2%
	Salinity (NaCl)	502	0-28.0	%	1%	±2%
	Specific Weight	503	1.000-1.220	-	0.001	±0.002
ORM 1SW	Refractive Index	504	0.0-99.9	%	0.1%	±0.2%
	Salinity Seawater	501	0-100	%	1%	±2%
	Chlorinity Seawater	502	0-37	%	1%	±2%
ORM 1AL	Specific Weight	503	1.000-1.070	-	0.001	±0.002
	Brix	504	0.0-50.0	%	0.1%	±0.2%
	Refractive Index	505	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1BR	Alcohol Meas.	501	0-72	%	1%	±1%
	Alcohol Vol.	502	0-80	%	1%	±1%
	Brix	503	0.0-50.0	%	0.1%	±0.2%
ORM 1BR	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Plato	501	0.0-20.0	°P	0.1	±0.3
	Brix	502	0.0-50.0	%	0.1%	±0.2%
ORM 1W	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Brix	501	0-150	°De	1	±2
	Vol%	502	0.0-22.0	%	0.1%	±0.2%
ORM 2W	KMW (Babo)	503	0.0-25.0	-	0.1	±0.2
	Brix	504	0.0-50.0	%	0.1%	±0.2%
	Vol%	502	0.0-22.0	%	0.1%	±0.2%
ORM 1CD	KMW (Babo)	503	0.0-25.0	-	0.1	±0.2
	Brix	504	0.0-50.0	%	0.1%	±0.2%
	Vol%	502	0.0-22.0	%	0.1%	±0.2%
ORM 2CD	Refractive Index	503	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	CorRee TDS 1	501	0.0-25.0	-	0.01	±0.20
	CorRee TDS 2	502	0.0-20.00	-	0.01	±0.20%
ORM 1UN	Refractive Index	503	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Brix	501	0.0-50.0	%	0.1%	±0.2%
	Specific Protein	502	0.0-20.0	g/100ml	0.1	±0.2
ORM 2UN	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Brix	501	0.0-50.0	%	0.1%	±0.2%
	Urea Cal	502	1.000-1.060	-	0.001	±0.002
ORM 1CA	Urea Cal	502	1.000-1.060	-	0.001	±0.002
	Urea Dng	503	0.0-50.0	%	0.1%	±0.2%
	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 2CA	Urea Cal	501	1.000-1.060	-	0.001	±0.002
	Brix	502	0.0-50.0	%	0.1%	±0.2%
	Refractive Index	503	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1C	Ethylenglycol (%)	501	0.0-100.0	%	0.1%	±0.5%
	Ethylenglycol (°C)	502	1.00-0.0	°C	0.1°C	±0.5°C
	Propylenglycol (%)	503	0.0-100.0	%	0.1%	±0.5%
ORM 2C	Propylenglycol (°C)	504	1.00-0.0	°C	0.1°C	±0.5°C
	Brix	505	0.0-99.0	%	0.1%	±0.2%