

Ring Illumination KERN OZB-IR

Professional illumination units guarantee outstanding, uniform and strong illumination

These illumination units are also available with UK mains plug. For more information on this, visit our online shop or give us a ring







OBB-A6102



OZB-A7101

Features

- Choose your favourite external illumination here to achieve maximum flexibility and greatest possible ease of use in stereo microscopy
- These professional illumination units provide a quality of light at a high, constant intensity at all times
- Regardless of whether your choice is spacesaving ring lights or cold light sources using optical fiber, our range is all you can wish for

OZB-A4572

- With the OZB-A7101 polarisation ring illumination unit, you also have an excellent component which has been specially optimised for observing shiny surfaces
- Naturally, these external illumination units also fit your standard stereo microscope
- Exception: The ring illumination units cannot be used in combination with the following ranges: OSE-1, OSF-4G, OZL-45R, OZC-5 and OZG-4

Model	Illuminance	Inner Ø	Colour temperature	Brightness adjustable	Illumination by segments	Polarising filter	
KERN		mm	K				
OZB-A4571	4W-LED	60	7000 - 11000	✓			
OZB-A4572	4W-LED	60	6500 – 7000	✓	✓		
OBB-A6102	4,5W-LED	63	approx. 7600	✓			
OZB-A7101	4,5W-LED	62	6500 – 7000	✓		✓	

✓ = Included with delivery

Goose Neck Illumination KERN OZB-IF







Application example

Features

 With the OZB-A4516 20 W-LED Goose Neck Illumination Unit with focusable light beam, you can adjust the illumination to suit your needs. Spot or scattered radiation means that you can achieve the very best illumination of your sample

✓
✓

✓ = Included with delivery

MICROSCOPES & REFRACTOMETERS 2024

KERN Pictograms





360° rotatable microscope head



Monocular MicroscopeFor the inspection with one eye



Binocular MicroscopeFor the inspection with both eyes



Trinocular MicroscopeFor the inspection with both eyes and the additional option for the connection of a camera



Abbe Condenser

With high numerical aperture for the concentration and the focusing of light



Halogen illumination For pictures bright and rich in contrast



LED illumination

Cold, energy-saving and especially long-life illumination



Incident illuminationFor non-transparent objects



Transmitting illuminationFor transparent objects



Fluorescence illumination For stereomicroscopes



Fluorescence illumination for compound microscopes

With 100 W mercury lamp and filter



Fluorescence illumination for compound microscopes

With 3 W LED illumination and filter



Phase contrast unit

For a higher contrast



Darkfield condenser/ unit

For a higher contrast due to indirect illumination



Polarising unit
To polarise the light



Infinity system Infinity corrected optical system



Zoom magnification For stereomicroscopes



·



Auto-focus

For automatic control of the focus level



Parallel optical system For stereomicroscopes, enables fatigue-proof working



Integrated scale

In the eyepiece



SD cardFor data storage



USB 2.0 interface For data transmission



USB 3.0 interface For data transmission



WIFI data interface:

For transmitting of the picture to a mobile display device



HDMI digital camera

For direct transmitting of the picture to a display device



PC software

To transfer the measurementsfrom the device to a PC.



Automatic temperature compesation

For measurements between 10 °C and 30 °C



Protection against dust and water splashes IPxx:

The type of protection is shown in the pictogram cf. DIN EN 60529:2000-09, IEC 60529:1989+A1:1999 +A2:2013



Battery operation

Ready for battery operation. The battery type is specified for each device.



Battery operation rechargeable

Prepared for a rechargeable battery operation



Plug-in power supply

230V/50Hz in standard version for EU.
On request GB, AUS or USA version.



Integrated power supply unit

Integrated in microscope. 230V/50Hz standard EU. More standards e.g. GB, AUS or USA on request.



Package shipment

The time required to manufacture the product internally is shown in days in the pictogram.



Pallet shipment

The time required to manufacture the product internally is shown in days in the pictogram.

Abbreviations

C-Mount Adapter for the connection of a

camera to a trinocular microscope

FPS Frames per second

H(S)WF High (Super) Wide Field (Eyepiece with high eye

point for wearers of glasses)

LWD Long Working Distance

N.A. Numerical Aperture

SLR camera Single-Lens Reflex camera

SWF Super Wide Field (Field number at least Ø 23 mm

for 10× eyepiece)

W.D. Working Distance

WF Wide Field (Field number up to Ø 22 mm

for 10× eyepiece)